The image features a background of a grid of green oval shapes that curves from the bottom left towards the top right. At the top left, the word "PODOSTRIPE" is written in white, bold, uppercase letters inside a green, slightly curved rectangular banner. At the bottom left, the word "PODOSTRIPE" is repeated in large, white, bold, uppercase letters, with a registered trademark symbol (®) to its upper right. Below this, five lines of text in white, bold, uppercase letters list the instruction manuals in different languages: German, English, Dutch, Italian, and Spanish. On the right side, there is a vertical column of four white rectangular buttons, each containing a language code in green, bold, uppercase letters: "DE", "EN", "NL", "IT", and "ES".

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Safety Instructions

Purpose:

The therapy with nail clips from 3TO GmbH serves to restore the natural growth of the nail. With the braces, pain can be relieved or the growth of ingrown and/or deformed nails can be corrected by lifting the ingrown nail. The brace is a disposable product and is delivered non-sterile. Nail correction braces are Class I non-invasive medical devices.

Indications:

Unguis incarnatus (ingrown nail), Unguis convolutus (rolled nail), Onychophosis (excessive hornification in the nail fold), Clavi in the nail fold (corn), after operations to prevent recurrence.

Contraindications:

The braces may only be used for the following clinical pictures after careful consideration of the benefit-risk ratio and, if necessary, in consultation with a doctor Onychomycosis, onycholysis, growth arrest, tumorous changes, subungual granulation, onychorrhexis, pAVK - peripheral arterial occlusive disease - (in individual cases), polyneuropathy / angiopathy (in individual cases), psoriasis (in individual cases), allergic reactions to ethyl cyanoacrylate.

Group of users:

To ensure the protection of the patient, the nail correction braces from 3TO GmbH may only be used by medically trained specialists (e.g. podiatrists or doctors). The medical expertise to assess the risks of the therapy as well as the ability to handle instruments are basic requirements for the application of nail correction braces. The liability for consequences of improper use shall be excluded herewith.

EN

Duration of the application:

The total duration of the therapy depends on the initial findings and the individual growth of the nail. The frequency of follow-up examinations with replacement of the tamponade fleece must be determined individually according to the medical history. If necessary, a daily check must be carried out. The braces usually have to be changed every 3 - 6 weeks.

Patient and user protection:

Clean and disinfect the nail to be treated. To protect the patient and user, the therapist must wear sterile, commercially available disposable gloves. For the therapist's own protection, he or she should wear a face mask and protective goggles during each treatment. Only use the braces as described in this manual and with original accessories.

Risks:

Risk of skin irritation, injury, inflammation and pain in the nail environment: The correct fitting and insertion of nail correction braces is an important prerequisite to avoid skin irritation, skin injury, pain, cornification and inflammation.

Risk of nail separation (onycholysis):

The correct adjustment of the pulling and levering forces of nail correction braces are important prerequisites to prevent the nail from loosening.

Risk of infection:

An infection of the nail plate, a pre-damage of the nail or a lack of hygiene (especially during repositioning) can lead to a nail mycosis. Patient education should include ongoing hygiene measures to prevent infection.

Danger from sticking:

The OrthoGlue adhesive bonds skin and eyelids together within seconds. Contains cyanoacrylate. Irritating to eyes, respiratory system and skin. Avoid contact with eyes and skin. In case of contact with eyes or internal Organs, please see your doctor immediately. Do not inhale vapours. May not come into contact with wounds - Cover wounds with oil or grease before use. Do not remove any skin adhesions by force but soak them in warm soapy water for a longer period of time or carefully loosen them with acetone.

The reactions usually occur quickly and are quickly noticed by the patient and user. Countermeasures must be taken immediately.

Special incidents:

Should any special incidents occur when using products of 3TO GmbH, remove the brace immediately from the patient. The affected toe must be disinfected and covered with sterile dressing. Consult a doctor immediately. Please contact 3TO GmbH immediately by phone +49 89 4520833-0 (outside opening hours by e-mail to info@3to-gmbh.de). You will then be sent a „Special incidents“ form. This is to be filled in and sent to 3TO GmbH as soon as possible.

Patient information:

Inform your patient about the causes and prevention possibilities of ingrown nails.

1 Remove non-adhesive parts from the nail surface. Smooth extremely rough nails or roughen extremely smooth nails **very lightly** with a file or burr.

2 Select the size so that the brace covers as much nail width as possible.

However, the expansion of the brace on the nail (approx. 0.7 - 2.7 mm) must be taken into account. Therefore, select the brace accordingly **shorter** than the nail width by the expansion length.

The brace should be placed as close as possible to the edge of the nail on the painful or more ingrown side.

3 Clean and degrease the brace and nail thoroughly with the alcohol swab and allow the alcohol to dry off well (**approx. 30 seconds**).

Health risk: danger of infection - see safety instructions


4 Open the adhesive and wipe the brush **lightly from both sides** of the bottle once at the bottle edge.

5 Wet the adhesive surface A of the pad with very little adhesive. Use the adhesive sparingly; if too much adhesive is used, the brace will not stick to the nail. Keep the edge and the thread of the bottle clean. Clean them if necessary, e.g. with a swab.

Health risk: Danger due to sticking - see safety instructions


EN

- 6** Position the adhesive surface A on the more ingrown side of the nail. Press the adhesive surface A against the nail with your finger for approx. **20 seconds** (do not move during this time!), then the adhesive must harden for **2 minutes** without further contact.


Health risk: Danger due to sticking - see safety instructions


- 7** Gently wipe the brush from both sides. Apply a thin layer of adhesive to the adhesive **surface B** of the brace. Make sure that no force is applied to the glued side A and that it is not pulled off.

Health risk: Danger due to sticking - see safety instructions


8 - 10 The brace is stretched by 0.7 to 2.7 mm (depending on the width of the brace, the nail thickness and the findings). Pull the brace in parallel position to the nail. Place the brace on the nail while stretched and press the adhesive surface B for **approx. 30 seconds**. Make sure that the brace is no longer displaced after it has been put on the nail. Afterwards the adhesive must harden for **2 minutes** without further touching.

Health risk: Danger due to sticking
see safety instructions

EN

11 Place the PODOSTRIPE pliers diagonally under the holding element and pinch off the pad.

To remove PODOSTRIPE, the edges of the brace are detached with the pliers. The brace can then be removed. Use a file or a burr to remove remains of the glue from the nail.

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